Produced By HuVeSearch BVBA/SPRL
Europalaan, 11 - 3900 PELT - BELGIUM

LEPEX – LEPEX Instruction Manual

The LEPEX extender is produced through the mixing and fine crushing of biochemical compounds (without water and being of pharmaceutical quality) and a minimal quantity of antibiotics in a controlled environment, with the help of special equipment and under special conditions of use. This process results in the formation of a high-quality, homogeneous, micronised powder which is easily soluble in water. After the preparation of the extender, it is essential to allow the mixture to stand for one hour to ensure complete hydration.

This LEPEX product is stored in polypropylene vials which are hermetically sealed with the help of a tight screw-on lid.

Each vial contains 62 grams of powder which corresponds to the preparation of one litre of the diluent.

Collecting, evaluating, diluting and storing the rabbit sperm

1. Preparing the extender

The extender can be prepared in a number of ways.

- a. Dilute the powder the day before collection in an adequate quantity of distilled water at room temperature and then store the mixture in the fridge. *If clumps* begin to form, these can be easily broken up using a spoon.* On the day of the collection, gently shake the extender and then store it at 32°C. The extender then needs to be kept in an oven and maintained at a temperature of 32°C.
- b. On the day of the collection, dilute the powder in the correct amount of water at 32°C. *If clumps* begin to form, these can be easily broken up using a spoon.* Leave to set for one hour before gently shaking. The extender then needs to be kept in an oven and maintained at a temperature of 32°C.

***Note: As a result of a physical reaction, small portions of micronised powder can clump. This is a natural occurrence and does not have any negative effect on the high quality of the extender. It is advisable to keep the extender in the fridge (4°C) to minimise the formation of clumps and also to respect the storage duration of the product.**

Any supplementary products (such as GnRH types - which are not delivered by HuVeSearch and are solely the responsibility of the user) being used should be added to the mixture shortly before the collection.

Note: Any surplus of the extender can be aliquoted into plastic tubes and placed in the freezer in order to be used in the weeks that follow. On the day of the collection, the tubes should be taken out of the freezer and placed in a water bath at 32°C. Alternatively, they can be placed in the fridge the day before collection to allow them to defrost slowly.

Produced By HuVeSearch BVBA/SPRL
Europalaan, 11 - 3900 PELT - BELGIUM

2. Preparing the artificial vagina.

The artificial vagina can be prepared in numerous ways.

- a. Prepare the artificial vaginas the day before collection. Fill with cold water and place them in an incubator at 45°C.
- b. Prepare the artificial vaginas on the day of collection. Fill them with water heated to approximately 50°C. It is essential to be sufficiently exact and to check the temperature with the help of a thermometer.

Bring the number of collection vials needed to the desired temperature. This can be done in two ways:

- a. Keep the collection vials in an incubator at 32°C with the extender (32 °C)
- b. Store the vials in an incubator or a water bath at 50 °C.

3. Collecting the sperm.

Shortly before carrying out the collection of sperm, place the collection vial in the collection system (artificial vagina). The time between set-up and collection needs to remain as short as possible in order to avoid the cooling down of the collection vial. Collect the sperm following normal procedures.

4. Macroscopic visual observation of the sperm

Observe the collected sperm in order to evaluate defects and quality.

If there is not enough sperm or if the sperm is too transparent, yellow, brown or red, the collected sperm cannot be used.

If any gel is present, this should be removed with the help of a pipette.

Good quality sperm has a milky, white appearance.

5. Initial dilution of sperm

After quickly observing the sperm, dilute the sperm in the vial using the extender. This can be done in a number of ways:

- a. Carefully add a volume of the extender equal to the amount of the collected sperm and shake carefully
- b. Carefully add a couple of millilitres of the extender to the sperm to dilute the sperm to a rough ratio of 1:10. (1 part sperm to 9 parts extender)

6. Microscopic evaluation of the sperm

For a microscopic evaluation of the sperm, the sperm should further be diluted. To do this, take a small quantity of pre-diluted sperm and dilute it again in the same extender in order to obtain a final dilution with a ratio of 1:200.

Produced By HuVeSearch BVBA/SPRL
Europalaan, 11 - 3900 PELT - BELGIUM

In the case of a first step 1:1 dilution (ref. point A above), it is necessary to dilute further to a ratio of 1:100.

In the case of a pre-dilution of 1:10 (ref. point B above), it is necessary to dilute the pre-diluted sperm further to a ratio of 1:20.

Take 10 micro-litres of this diluted sperm and place onto a microscopic slide (preheat to 37°C) and place a 22x22mm cover-slide over it.

A counting chamber can also be used.

Only use ejaculates which have a total mobility of above 75%.

7. Timing the mixing of different ejaculates.

The desired end result is the mixing of different ejaculates in order to be able to distribute high-quality, commercial quantities. This can be achieved in a number of different ways:

- a. Mixing sperm directly after ejaculation
- b. Mixing pre-diluted sperm
- c. Mixing post-diluted sperm which has the final desired concentration.

Each method has both advantages and disadvantages.

Method A appears to be the simplest but prevents the microscopic testing of different samples and forces the use of a fixed-ratio dilution.

Method B also prevents microscopic testing but works on the assumption that it is better to allow the cells to stabilise in a smaller quantity of extender.

Method C has the advantage of being able to optimise each ejaculate individually after microscopic observation.

In all circumstances, it is essential to proceed in the quickest amount of time possible to the dilution steps, including the last one, to guarantee the quality of sperm.

8. Final dilution of the sperm

The final dilution of the pre-diluted sperm occurs in order to optimise the production on the basis of the final concentration of sperm cells. The final concentration can be influenced by different factors. A target concentration of 20 million cells per millilitre acts as a good reference value. For precise calculations, the only solution is the use of a counting chamber.

9. Storage of diluted sperm

Once it has reached its target concentration, the diluted sperm should be stored at an ideal temperature of 17°C. In order to maintain fertility, it is essential that the time between the collection of the sperm and its storage at 17°C is reduced as much as possible.

Produced By HuVeSearch BVBA/SPRL
Europalaan, 11 - 3900 PELT - BELGIUM

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