



**SERGA SL**  
RABBIT BREEDING, SELECTION  
AND ARTIFICIAL  
INSEMINATION CENTER

Distributed by HuVeSearch BV  
Europalaan 11 B8  
3900 Pelt - Belgium



# DATASHEET

## ARTIFICIAL VAGINA FOR RABBIT SEMEN COLLECTION



Flexible and solid devices for rabbit semen extraction.

### CHARACTERISTICS

- **Efficient Collection:** Saves time during extraction
- **Greater Safety:** Eliminates the need to warm the semen with water and the associated risks.
- **Increased Ejaculate Volume:** Optimized anatomical design.
- **Improved Rabbit Welfare:** Prioritizes comfort to reduce stress.

### THECNICAL SPECIFICATIONS

|                     |                           |
|---------------------|---------------------------|
| Material:           | biomedical-grade silicone |
| Shape:              | ergonomic design          |
| Height:             | 80 mm - 3 1/8"            |
| Widht:              | 43 mm - 1 3/4"            |
| Depth:              | 43 mm - 1 3/4"            |
| Weight:             | 111 g                     |
| Thermic resistance: | 110°C - 231,8°F           |
| Paint:              | pigment for silicones     |
| Properties:         | hypoallergenic            |
|                     | high tear resistance      |

### USE

- Heat to a maximum of 80°C in an isolated and sterile environment.
- Operating temperature: 45°C at the time of semen collection.
- Collect using a conical PP tube with a diameter of 16 mm.

### CLEANING

Use of mild soap.  
Disinfection in a sterilizing or drying oven at 100°C for 30 minutes.

### COLORS

**Blue:** Sensitive Model,  
extra elasticity/stretch.

**Pink:** Comfort Model,  
high elasticity/stretch.

Please visit our website: [www.lepex-novum.org](http://www.lepex-novum.org)